



Recombinant Danio rerio Transcriptional adapter 3 (tada3)

Product Code	CSB-MP742457DIL
Abbreviation	tada3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SY21
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSELKDCPPL KYYDFKPVDH VKVCPRYTAV LSRSEDDGIG IEELDTLQLE LETLLSSASR RLRALEEQRQ ILTDWQDKKG DKRFLKLEKD PDLAASSRHS KPKKQKLDVK GSHGPGPGPG RPKSKNIQTK VQDFEFEVDP QDIPRNPKNL APNRFWASVE PYCADITNEE IRVLEELLKA PDDEAEYYKI PTLGKHYSQR WAQEDLLEER REGARANDKK KSMMGPLSEL DAKDVDALLK KSESQHEPPE DGCPFGPLSQ RLLQALVEEN IISPMEDSPI PDIPGKDDGA GTSPRSQGKA FSVPHTRSLE ARIREELVSQ GLLDSDERQG VGGESEDEVL AELQKRQAEK KALTAHNRSR KLELLKLARE EMRKQELRQR VRVADNEVME AFRRIMAARQ KKRTPTKKEK DQAWKALKER ESILKLLDG
Source	Mammalian cell
Target Names	tada3
Protein Names	Recommended name: Transcriptional adapter 3 Alternative name(s): ADA3 homolog Transcriptional adapter 3-like Short name= ADA3-like protein
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.