



Recombinant Danio rerio Tumor necrosis factor alpha-induced protein 8-like protein 1 (tnfaip8l1)

Product Code	CSB-BP742464DIL
Abbreviation	tnfaip8l1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SZE8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDSFSTKNLA LQAQKKLMSK MATKTVANLF IDDTSSSEVLD ELYRVTKEYT RNRKEAQKII KNLIKMVVKL GVLYRNGQFN NEELALVERF RKKVHTLAMT AVSFYQIDFT FDRRVMSNLL NDCRELLHQA INRHLLAKSH ARINHVFNFH ADCDFLATLY GPSEVYRGHL QKICEGVNKM LDEGNL
Source	Baculovirus
Target Names	tnfaip8l1
Protein Names	Recommended name: Tumor necrosis factor alpha-induced protein 8-like protein 1 Short name= TNF alpha-induced protein 8-like protein 1 Short name= TNFAIP8-like protein 1
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.