



Recombinant Danio rerio Dynactin subunit 2 (dctn2)

| | |
|--------------------------|--|
| Product Code | CSB-EP742480DIL-B |
| Abbreviation | dctn2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q7T3H1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MADPKYANLP GIASNEPDVY ETSIDLPEDDQ AQFESELEEL CSDSVERIVV NPNAAYDKFK DKHVSASLD FSDRISKNNR VGYESGDYEL LAEGCGVKET PQQKYQRLVN EIHLCQDVE KIQTSTKESG AEERLTPVAL AQQAAQLKQQ LVSAHLDSL L GPDHINLTD PDGALAKRLL TQLEVARGVR SGAGADGKTA APKGPDGVL YELHSRPEQE KFTESAKVAE LERRLALET AVGSGSDKRG PLSSGVQGSS LTETLELLQA RVSALDAATL DQVEARLQSV LGKMNEIAKH KATMEDAETQ SKVSQLYDVV QKWDAMATSL PQVVRRLTAV RELHEQAMQF GQLLTHLDTT QQMINNSLKD NSTLLTQVQQ TMKENLLAVE ENFSALDQRM KKLNK |
| Source | E.coli |
| Target Names | dctn2 |
| Protein Names | Recommended name: Dynactin subunit 2 |
| Expression Region | 1-405 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |