



# Recombinant Mouse Zinc finger protein ZIC 5 (Zic5)

<b>Product Code</b>	CSB-EP742510MO
<b>Abbreviation</b>	Zic5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7TQ40
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMEPPLSKRN PPALRLADLA TAQAQQLQNM TGFPVLVGPP AHSQRRVAVM HLHPRDLGTD PGVASTALGP EHMAQASGQG PCPPSQGLPG LFQVPAPAAR SVASGTHPGA RTHPDGGGSS GAQASAPPPP APPLPPSQSS SPPPPPPPPP ALSGYTATNS GGGSSSGKGH SRDFVLRRLD SATAPAAAMH GAPLGGEQRS GSSSPQHPTP PPHPAGMFIS ASGTYAGRDG GGSALFPALH DSPGAPGGHP LNGQMRLGLA AAAAAAAELY GRAEPPFAPR SGDAHYGAVA AAAAAALHGY GAVNLNLNLA AAAAAAAAAG PGPHLQHHAP PPAPPAPAP HPHHPHLPGA AGAFLRYMRQ PIKRELICKW LDPEELAGPP ASADSGVKPC SKTFGTMHEL VNHVTVEHVG GPEQSSHVCF WEDCPREGKP FKAKYKLINH IRVHTGEKPF PCFPFGCGKV FARSENKIH KRTHTGEKPF KCEFDGCDRK FANSSDRKKH SHVHTSDKPY YCKIRGCDKS YTHPSSLRKH MKIHCKSPPP SPGALGYSSV GTPVGDPLSP VLDPTRSRSS TLSPQVTNLN EWYVCQASGA PSHLHTPSSN GTTSESEDEE MYGNPEVMRT IH
<b>Source</b>	E.coli
<b>Target Names</b>	Zic5
<b>Protein Names</b>	Recommended name: Zinc finger protein ZIC 5 Alternative name(s): Zinc finger protein of the cerebellum 5
<b>Expression Region</b>	1-622
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the ZIC family of C2H2-type zinc finger proteins. Members of this family are important during development, and have been associated X-linked visceral heterotaxy and holoprosencephaly type 5. This gene is closely linked to a gene encoding zinc finger protein of the cerebellum 2, a related family member on chromosome 13.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.