



Recombinant UPF0232 protein Mb0004 (Mb0004)

Product Code	CSB-EP742566MVH-B
Abbreviation	Mb0004
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7U313
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
Purity	>85% (SDS-PAGE)
Sequence	MTGSVDRPDQ NRGRLMKSP GLDLVRRRLD EAAAAARARG QDAGRGRVAS VASGRVAGRR RSWGPGPDI RDPQPLGKAA RELAKKRGWS VRVAEGMVLG QWSAVVGHQI AEHARPTALN DGVLSVIAES TAWATQLRIM QAQLLAKIAA AVGNDDVRSR KITGPAAPSW RKGPRHIAGR GPRDTYG
Source	E.coli
Target Names	BQ2027_MB0004
Protein Names	Recommended name: UPF0232 protein Mb0004
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.