



Recombinant tRNA pseudouridine synthase A (truA)

| | |
|--------------------------|---|
| Product Code | CSB-MP742571BIK |
| Abbreviation | truA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q7U376 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bordetella parapertussis (strain 12822 / ATCC BAA-587 / NCTC 13253) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSRIALGLAY DGS AWQGWQT QPHGVTVQDQ VEAALASFAG GGGPVATVCA GRTDTGVHAA MQVIHLDTDL QRRDESWVRG VNAFLPPSIV VQWARPVSEA FHARFSARSR TYVYLLWRGR VRPALWAGRA GWAFQPLDVP AMRAAARTLL GEHDFSSFERS SQCQARHPVR TLHRLDIDER GAFLVFTLRA NAFLHHMVRN LIGALLQVGQ GRESVAWMDA LLRARDRRLG APTFMPDGLY LSAIEYPAEF GFDEL DGGTM LLS PFTGALG |
| Source | Mammalian cell |
| Target Names | truA |
| Protein Names | Recommended name: tRNA pseudouridine synthase A EC= 5.4.99.12 Alternative name(s): tRNA pseudouridine(38-40) synthase tRNA pseudouridylate synthase I tRNA-uridine isomerase I |
| Expression Region | 1-270 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |