



Recombinant Staphylococcus aureus Uncharacterized leukocidin-like protein 2 (SAR2108)

Product Code	CSB-EP743339SKV
Abbreviation	SAR2108
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GF49
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Purity	>85% (SDS-PAGE)
Sequence	DSQDQTKK EHVDKAQQKE KRNVDKDKN TPGPDDIGKN GKVTKRTVSE YDKETNILQN LQFDFIDDP YDKNVLLVKK QGSIHNLKF ESHRNETNAS WLKYPSEYHV DFQVQRNPKT EILDQLPKNK ISTAKVDSTF SYSLGGKFDS TKGIGRTSSN SYSKSISYNQ QNYDTIASGK NNNRHVHWSV VANDLKYGNE IKNRNDEFLF YRNTRLSTVE NPELSFASKY RYPALVRSGF NPEFLTYISN EKSNEKTRFE VTYTRNQDIL KNKPGIHYGQ PILEQNKDGQ RFIVVYEVDW KNKTVKVVEK YSDQNKPYKE G
Source	E.coli
Target Names	SAR2108
Protein Names	Recommended name: Uncharacterized leukocidin-like protein 2
Expression Region	33-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.