



# Recombinant *Oryza sativa* subsp. japonica Putative protein phosphatase 2C 24 (Os02g0633900, LOC\_Os02g42270)

<b>Product Code</b>	CSB-YP743478OFG
<b>Abbreviation</b>	Os02g0633900, LOC_Os02g42270
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6H7J3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEALPQIRQT LSEIDRRIPD ALRVAMGLRL RPTAGAALEE VTRIAASCLP RPCPEGGDDP MECDEAAPAR ALRMEAASCF LPDHDEDTHF VRPEAGVVAL ADGVGGYRAP GVDAAAFARA LMYNAFEMVV ATTPGGAGGI CPYALLGWAY EQAVSARTQG ASTAVILSLA GATLKYAYIG DSAFAVFRDG KLFFRSEAQV HSFNYPFQLS VKNGNSV TSA ARGGVEVKEG DVVVAGTDGL FDNVTSEELQ RIVAMGRALG LSPKQTADV AGFAYEASTT MGRDTPFSLE SRKKQGTIFR RGKRDDITVV VAYIV
<b>Source</b>	Yeast
<b>Target Names</b>	Os02g0633900
<b>Protein Names</b>	Recommended name: Putative protein phosphatase 2C 24 Short name= OsPP2C24 EC= 3.1.3.16
<b>Expression Region</b>	1-315
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.