



Recombinant Human Transport and Golgi organization protein 2 homolog (TANGO2)

Product Code	CSB-BP743560HU
Abbreviation	C22orf25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ICL3
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MCIFFKFDP RPVSKNAYRL ILAANRDEFY SRPSKLADFW GNNNEILSGL DMEEGKEGGT WLGISTRGKL AALTNYLQPQ LDWQARGRGE LVTHFLTTDV DSLSYLKKVS MEGHLYNGFN LIAADLSTAK GDVICYYGNR GEPDPIVLTP GTYGLSNALL ETPWRKLCFG KQLFLEAVER SQALPKDVLI ASLLDVLNNE EAQLPDPAIE DQGGEYVQPM LSKYAAVCVR CPGYGTRTNT IILVDADGHV TFTERSMMDK DLSHWETRTY EFTLQS
Source	Baculovirus
Target Names	TANGO2
Protein Names	Recommended name: Uncharacterized protein C22orf25
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.