



Recombinant *Oryza sativa* subsp. japonica Auxin-responsive protein IAA9 (IAA9)

Product Code	CSB-MP743618OFG
Abbreviation	IAA9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6K846
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MELELGLAPP NSGHLVVDEL SSSSSSGGGS GSAPVSASSA GKRGFREAFAQ ETLLLFDDGS CCNTSDDDCR RRKKTVVGWP PVSSARRACG GANYVKVKKE GDAIGRKVDL ALHSSYDELA ATLARMFPTN DHQGEKKMAN DDHGDAAGPV VTYEDGDGDW MLVGDVPWDD FARSVKRLKI LG
Source	Mammalian cell
Target Names	IAA9
Protein Names	Recommended name: Auxin-responsive protein IAA9 Alternative name(s): Indoleacetic acid-induced protein 9
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.