



# Recombinant *Oryza sativa* subsp. japonica Probable E3 ubiquitin-protein ligase XBOS32 (XBOS32)

<b>Product Code</b>	CSB-YP743621OFG
<b>Abbreviation</b>	XBOS32
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6KAE5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGFLSLVGNS FGCSASGERL VSAARDGDLQ EARALLEYNP RLARYSTFGG RNSPLHYAAA QGHHEIVSLL LESGVEINLR NYRGQTALMQ ACQYGHWEVV QTLMLFNANV HRTDYLNNGS ALHFAALHGH ARCLRLVLAD YVPSMPNFWN SMKDSLSEEG PSADLDEDGL FKMVNQKADG GLTPLHMAAL NGHVECVQLL LDLGASVIEA TIEDGTTIDL IGAGSTPLHY AACGGNAVCC QLLIARGASL SAQNASGWTP LMVARSWHRN SLEEILSKEP ESRI RTVPSP YLCLPLMSIM SIAREFGWRY LNQSPVCIDP CAVCLEGSCS VAAEGCKHEF CTRCALYLCS TSYTSVSPAG AIPCPLCRHP IIAFTALPGT SPIRELPRNS LLSLFCCTTCP AVNSDSTPSI ASHLRYRTEFQ CARMPPMGSS SFRSLSCQRL PAMKLNPSFC MGAMDTNPCL IRCSRFGPSF RRSASQGESS RRAWPLTFDP IAATGS
<b>Source</b>	Yeast
<b>Target Names</b>	XBOS32
<b>Protein Names</b>	Recommended name: Probable E3 ubiquitin-protein ligase XBOS32 EC= 6.3.2.- Alternative name(s): Ankyrin repeat domain and RING finger-containing protein XBOS32 XB3 protein homolog 2
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.