



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 39 (Os05g0576300, LOC_Os05g50080)

Product Code	CSB-EP743686OFG
Abbreviation	Os05g0576300, LOC_Os05g50080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6L5G1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MQEALLPPAH PGRFYSDFGP KPFSGDQRL SSPNLLTNGG DLFYGCYSPF SPTRVLSPPP PRRAASFSHC SSSSDSVVDD GDGAGAAAAT EHLHLAHLA LQYQEMANRF ELCLSHLADA ADEAAALRQE NAELRVANND LACRIAKFGG RQSSAIALAG DLRRRLRPKE QTVPALPPPP QSPPAALMNP VAVPEKQAVL PKSISIRSTG YQKLNQGGKH RVSKPVNVGS QRVFVIGIDGA EGGEHKVGVK KEEPPMGGLE FEVYNQGMFK TELCNKWEET GACPYGDQCQ FAHGVAELRP VIRHPRYKTQ VCRMVLAGGV CPYGHRCHFR HSITPADRFS FGH
Source	E.coli
Target Names	Os05g0576300
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 39 Short name= OsC3H39
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.