



# Recombinant Danio rerio AP-2 complex subunit mu-A (ap2m1a)

<b>Product Code</b>	CSB-MP744006DIL
<b>Abbreviation</b>	ap2m1a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NWK2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIGGLFIYNH KGEVLISRVY RDDIGRNAVD AFRVNVIHAR QQVRSPVTNI ARTSFFHVKR SNIWLAAVTK QNVNAAMVFE FLYKMCDVMT AYFGKISEEN IKNNFVLIYE LLDEILDFGY PQNSETGALK TFITQQGIKG QHQTKEEQSQ ITSQVTGQIG WRREGIKYRR NELFLDVLES VNLLMSPQGQ VLSAHVSGRV VMKSYLSGMP ECKFGMNDKI VIDKQKGKGGT TDDAGKSGKQ SIAIDDCTFH QCVRLSKFDS ERSISFIPPD GEYELMRYRT TKDIILPFRV IPLVREVGRT KLEVKVIKS NFKPSLLAQK IEVRIPTPLN TSGVQVICMK GKAKYKASEN AIVWKIKRMA GMKESQISAE IELLPTNDKK KWARPPISMN FEVPFAPSGL KVRYLKVFES KLNYSDDHVI KWVRYIGRSG IYETRC
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ap2m1a
<b>Protein Names</b>	Recommended name: AP-2 complex subunit mu-A Alternative name(s): AP-2 mu-A chain Clathrin assembly protein complex 2 medium chain A Clathrin coat assembly protein AP50-A Clathrin coat-associated protein AP50-A Mu2-adaptin-A Pla
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.