



Recombinant Human Putative deoxyribonuclease TATDN1 (TATDN1)

Product Code	CSB-EP744031HU-B
Abbreviation	TATDN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P1N9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSRFKFIDIG INLTDPMFRG IYRGVQKHQD DLQDVIGRAV EIGVKKFMIT GGNLQDSKDA LHLAQTNGMF FSTVGCHPTR CGEFEKNNPD LYLKELLNLA ENNKGVVAI GECGLDFDRL QFCPKDTQLK YFEKQFELSE QTKLPMFLHC RNSHAEFLDI MKRNRDRCVG GVVHSFDGTK EAAAALIDL LYIGFNGCSL KTEANLEVLK SIPSEKLMIE TDAPWCGVKS THAGSKYIRT AFPTKKKWES GHCLKDRNEP CHIIQILEIM SAVRDEDPLE LANTLYNNTI KVFFPGI
Source	E.coli
Target Names	TATDN1
Protein Names	Recommended name: Putative deoxyribonuclease TATDN1 EC= 3.1.21.- Alternative name(s): Hepatocarcinoma high expression protein
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.