



Recombinant Human RNA demethylase ALKBH5 (ALKBH5)

Product Code	CSB-MP744062HU
Abbreviation	ALKBH5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P6C2
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAASGYTDL REKLKSMTSR DNYKAGSREA AAAAAA AVAA AAAAAAAAEP YPVSGAKRKY QEDSDPERSD YEEQQLQKEE EARKVKSGIR QMRLFSQDEC AKIEARIDEV VSRAEKGLYN EHTVDRAPLR NKYFFGEGEY YGAQLQKRG GQERLYPPGD VDEIPEWVHQ LVIQKLVEHR VIPEGFVNSA VINDYQPGGC IVSHVDPIHI FERPIVSVSF FSDSALCFGC KFQFKPIRVS EPLVSLPVRR GSVTVLSGYA ADEITHCIRP QDIKERRAVI ILRKTRLDAP RLETKSLSSS VLPPSYASDR LSGNNRDPAL KPKRSHRKAD PDAHRPRIL EMDKEENRRS VLLPTHRRRG SFSSSENYWRK SYESSEDCSE AAGSPARKVK MRRH
Source	Mammalian cell
Target Names	ALKBH5
Protein Names	Recommended name: RNA demethylase ALKBH5 EC= 1.14.11.-Alternative name(s): Alkylated DNA repair protein alkB homolog 5
Expression Region	2-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.