



Recombinant Rat Plasma protease C1 inhibitor (Serping1)

Product Code	CSB-BP744067RA
Abbreviation	Serping1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P734
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	DSEVTSHS SQDPLVVQEG SRDSVPERDG SRSPIEHTGQ SSTWPTTSGS TKISNDTMDQ VANESFIQHV QPAAQLPEDS PSQSPVNSSS PPSTASAPPT QAPTEPLCPE PLAWCSDSDR DSSEATLSEA LTDFSVKLYH AFSATKKAET NMAFSPFSIA SLLTQVLLGA GDSTKSNLED ILSYPKDFAC VHQTLKAFSS KGVTSVSQIF HSPDLAIRDT YVNASLSLYG SSPRVLGPDG DANLKLINTW VAENTNHKIN ELLDSLPSDT RLVLLNAVYL SAKWKKTFEQ KKMMASFLYK NSMIKVPMLS SKKYPLALFN DQTLKAKVGQ LQLSHNLSFV IMVPQSPTHQ LEDMEKALNP TVFKAILKKL ELSKFQPTYV MMPRIKVKSS QDMLSIMEKL EFFDFTYDLN LCGLTEDPDL QVSSMKHETV LETTETGVEA AAASISVAR NLLIFEVQQP FLFLLWDQRH KFPVFMGRVY DPRA
Source	Baculovirus
Target Names	Serping1
Protein Names	Recommended name: Plasma protease C1 inhibitor Short name= C1 Inh Short name= C1Inh Alternative name(s): C1 esterase inhibitor C1-inhibiting factor Serpin G1
Expression Region	23-504
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a highly glycosylated plasma protein involved in the regulation of the complement cascade. Its protein inhibits activated C1r and C1s of the first complement component and thus regulates complement activation. Deficiency of this protein is associated with hereditary angioneurotic oedema (HANE). Alternative splicing results in multiple transcript variants encoding the same isoform.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.