



# Recombinant Human Peptidyl-prolyl cis-trans isomerase CWC27 homolog (CWC27)

<b>Product Code</b>	CSB-YP744267HU
<b>Abbreviation</b>	CWC27
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6UX04
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           SNIIYIQEPP TNGKVLLKTT AGDIDIELWS KEAPKACRNF IQLCLEAYYD            NTIFHRVVPV FIVQGGDPTG TGSGGESIYG APFKDEFHSR LRFNRRGLVA            MANAGSHDNG SQFFFTLGRA DELNNKHTIF GKVTGDTVYN MLRLSEVDID            DDERPHNPHK IKSCEVLFNP FDDIIPREIK RLKKEKPEEE VKKLKPKGTK            NFSLLSFGEE AEEEEEEVNR VSQSMKGKSK SSHDLLKDDP HLSSVPVVES            EKGDAPDLVD DGEDESAEHD EYIDGDEKNL MRERIAKKLK KDTSANVKS            A GEGEVEKKS SVSRSEELRKEA RQLKRELLAA KQKKVENAAK QAEKRSEEEE            APPDGAVAEY RREKQKYEAL RKQQSKKGTS REDQTLALLN QFKSKLTQAI            AETPENDIPE TEVEDDEGWM SHVLQFEDKS RKVKDASMQR SDTFEYDPR            NPVNRKRRREE SKKLMREKKE RR         </p>
<b>Source</b>	Yeast
<b>Target Names</b>	CWC27
<b>Protein Names</b>	Recommended name: Peptidyl-prolyl cis-trans isomerase CWC27 homolog Short name= PPlase CWC27 EC= 5.2.1.8 Alternative name(s): Antigen NY-CO-10 Serologically defined colon cancer antigen 10
<b>Expression Region</b>	2-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.