



Recombinant Human Leucine-rich repeat-containing protein 31 (LRRC31)

Product Code	CSB-EP744274HU-B
Abbreviation	LRRC31
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UY01
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSQTRKKTSS EGETKPQTST VNKFLRGSNA ESRKEDNDLK TSDSQPSDWI QKTATSETAK PLSSEMEWRS SMEKNEHFLQ KLGKKAVNKC LDLNNCGLTT ADMKEMVALL PFLPDLEELD ISWNGFVGGT LLSITQQMHL VSKLKILRLG SCRLTTDDVQ ALGEAFEMIP ELEELNLSWN SKVGGNPLI LQKFQKGSKI QMIELVDCSL TSEDGTFLGQ LLPMLQSLEV LDLSINRDIV GSLNSIAQGL KSTSNLKV LK LHSCGLSQKS VKILDAAFY LGELRKL DLS CNKDLGGGFE DSPAQLVMLK HLQVLDLHQC SLTADDVMSL TQVIPLLSNL QELDLSANKK MGSSSENLLS RLRFLPALKS LVINNCALES ETFTALAEAS VHLSALEVFN LSWNKCVGGN LKLLLETLKL SMSLQVLRSL SCSLVTEDVA LLASVIQTGH LAKLQKLDLS YNDSICDAGW TMFCQNVRF LKELIELDISL RPSNFRDCGQ WFRHLLYAVT KLPQITEIGM KRWILPASQE EELECFDQDK KRSIHFDHGG FQ
Source	E.coli
Target Names	LRRC31
Protein Names	Recommended name: Leucine-rich repeat-containing protein 31
Expression Region	1-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.