



# Recombinant *Gossypium hirsutum* Tubulin alpha-2 chain

<b>Product Code</b>	CSB-EP744282GHB-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6VAG0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Gossypium hirsutum</i> (Upland cotton) ( <i>Gossypium mexicanum</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRECISIHIG QAGIQVGNAC WELYCLEHGI QPDGQMPSDK TVGGGDDAFN TFFSETGAGK HVPRAVFDL EPTVIDEVRT GAYRQLFHPE QLISGKEDAA NNFARGHYTI GKEIVDLCLD RIRKLADNCT GLQGFLVFNA VGGGTGSGLG SLLLERLSVD YGKKSCLGFT VYPSPQVSTS VVEPYNSVLS THSLEHTDV AVLLDNEAIY DICRRSLDIE RPTYTNLNL VSQVISSLTA SLRFDGALNV DVTEFQTNLV PYPRIHFMLS SYAPVISA EK AYHEQLSVAE ITNSAFEPSS MMAKCDPRHG KYMACCLMYR GDVVPKDVNA AVASIKTKRT IQFVDWCPTG FKCGINYQPP TVVPGGDLAK VQRAVCMISN STSVAEVFSR IDHKFDLMYA KRAFIHWYVG EGMEEGEFSE AREDLAALEK DYEEVGAESA EGDEGEDDEY
<b>Source</b>	E.coli
<b>Protein Names</b>	Recommended name: Tubulin alpha-2 chain Alternative name(s): Alpha-2-tubulin
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.