



Recombinant Human Developmental pluripotency-associated protein 3 (DPPA3)

Product Code	CSB-BP744290HU
Abbreviation	DPPA3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6W0C5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDPSQFNPTY IPGSPQMLTE ENSRDDSGAS QISSETLIKN LSNLTINASS ESVSPLSEAL LRRESVGA AV LREIEDEWLY SRRGVRTLLS VQREKMARLR YMLLGGVIRTH ERRPTNKEPK GVKKESRPFK CPCSFVCSNG WDPSENARIG NQDTKPLQP
Source	Baculovirus
Target Names	DPPA3
Protein Names	Recommended name: Developmental pluripotency-associated protein 3 Alternative name(s): Stella-related protein
Expression Region	1-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein that in mice may function as a maternal factor during the preimplantation stage of development. In mice, this gene may play a role in transcriptional repression, cell division, and maintenance of cell pluripotentiality. In humans, related intronless loci are located on chromosomes 14 and X.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.