



Recombinant Zea mays Fructokinase-2 (FRK2)

Product Code	CSB-YP744332ZAX
Abbreviation	FRK2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q6XZ78
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	≥85% (SDS-PAGE)
Sequence	MAPLGDGGAA AAAASNLLVV SFGEMLIDFV PDVAGLSLAE SGGFVKAPGG APANVACAIA KLGSSAFVG KFGDDEFHGM LVNILKQNNV NAEGCLFDKH ARTALAFVTL KHDGEREFMF YRNPSADMLL TEAELDGLLV RRARVFHYGS ISLISEPCRS AHMAAMRAAK AAGVLCSYDP NVRLPLWPSP DAAREGILSI WKEADFIKVS DDEVAFLTRG DANDEKNVLS LWFDGLKLLV VTDGDKGCRY FTKDFKGSVP GFKVDTVDTT GAGDAFVGSL LVNVAKDDSI FHNEEKLEA LKFSNACGAI CTTKKAIPA LPTVATAQDL IAKAN
Source	Yeast
Target Names	FRK2
Protein Names	Recommended name: Fructokinase-2 EC= 2.7.1.4 Alternative name(s): ZmFRK2
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>