



Recombinant Rat Myeloid differentiation primary response protein MyD88 (Myd88)

Product Code	CSB-MP744333RA
Abbreviation	Myd88
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Y1S1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MSAGGPRVGS VSVDSYLFSL PLVALNVGVR RRLSLFLNPR TTAAADWTSL AEEMGFEYLE IREFETRPDP TRSLLDWQGG RSGSSVGRLL ELLALLDRED ILYELKDRIE EDCQKYIRNQ QKQESEKPLQ VARVESSVPQ TKELGGITTL DDPLGQTPEL FDAFICYCPS DIEFVQEMIR QLEQTDYRLK LCVSDRDVLP GTCVWSIASE LIEKRCRRMV VVVSDDYLQS KECDFQTKFA LSLSPGVQQK RLIPIKYKAM KKDFPSILRF ITICDYTNPC TKSFWWTRLA KALSLP
Source	Mammalian cell
Target Names	Myd88
Protein Names	Recommended name: Myeloid differentiation primary response protein MyD88
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.