



Recombinant *Oryza sativa* subsp. japonica Molybdopterin synthase sulfur carrier subunit (Os02g0558300, LOC_Os02g35200)

Product Code	CSB-MP744362OFG
Abbreviation	Os02g0558300, LOC_Os02g35200
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YVX4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MALDPKANHA AAAAASADNP TAAAKAKVK VKVLFFARAR DLTGVTEAPV EVPAGSTAGD CLARVLAAFP RLEEIRSMV LALNEEYAPE DAAVGDGDEL AIIPPISGG
Source	Mammalian cell
Target Names	Os02g0558300
Protein Names	Recommended name: Molybdopterin synthase sulfur carrier subunit Alternative name(s): Molybdenum cofactor synthesis protein 2 small subunit Molybdenum cofactor synthesis protein 2A Short name= MOCS2A Sulfur carrier protein MOCS2A
Expression Region	1-109
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.