



Recombinant Human WD repeat-containing protein 88 (WDR88)

Product Code	CSB-EP744395HU
Abbreviation	WDR88
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZMY6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASPPRCSPT AHDRECKLPP PSAPASEYCP GKLSWGTMAR ALGRFKLSIP HTHLLATLDP LALDREPPPH LLPEKHQVPE KLIWGDQDPL SKIPFKILSG HEHAVSTCHF CVDDTKLLSG SYDCTVKLWD PVDGSVVRDF EHRPKAPVVE CSITGDSSRV IAASYDKTVR AWDLETGKLL WKVRYDTFIV SCKFSPDGKY VVSGFDVDHG ICIMDAENIT TVSVIKDHHT RSITSCCFDP DSQRVASVSL DRCIKIWDVT SQATLLTITK AHSNAISNCC FTFSGHFLCT SSWDKNLKIW NVHTGEFRNC GACVTLMQGH EGSVSSCHFA RDSSFLISGG FDRTVAIWDV AEGYRKLSLK GHNDWVMDVA ISNNKKWILS ASKDRMTMLW NIEEIDEIPL VIKYKKAVGL KLKQCERCDR PFSIFKSDTS SEMFTQCVFC RIDTRGLPAD TSSSSSSSER ENSPPPRGSK DD
Source	E.coli
Target Names	WDR88
Protein Names	Recommended name: WD repeat-containing protein 88 Alternative name(s): PQQ repeat and WD repeat-containing protein
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.