



# Recombinant Human WD repeat-containing protein 88 (WDR88)

<b>Product Code</b>	CSB-MP744395HU
<b>Abbreviation</b>	WDR88
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ZMY6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASPPRCSPT AHDRECKLPP PSAPASEYCP GKLSWGTMAR ALGRFKLSIP HTHLLATLDP LALDREPPPH LLPEKHQVPE KLIWGDQDPL SKIPFKILSG HEHAVSTCHF CVDDTKLLSG SYDCTVKLWD PVDGSVVRDF EHRPKAPVVE CSITGDSSRV IAASYDKTVR AWDLETGKLL WKVRYDTFIV SCKFSPDGKY VVGFDVDHG ICIMDAENIT TVSVIKDHHT RSITSCCFDP DSQRVASVSL DRCIKIWDVT SQATLLTITK AHSNAISNCC FTFSGHFLCT SSWDKNLKIW NVHTGEFRNC GACVTLMQGH EGSVSSCHFA RDSSFLISGG FDRTVAIWDV AEGYRKLKSLK GHNDWVMDVA ISNNKKWILS ASKDRMTMLW NIEEIDEIPL VIKYKKAVGL KLKQCERCDCR PFSIFKSDTS SEMFTQCVFC RIDTRGLPAD TSSSSSSSER ENSPPPRGSK DD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WDR88
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 88 Alternative name(s): PQQ repeat and WD repeat-containing protein
<b>Expression Region</b>	1-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.