



# Recombinant Human Coiled-coil domain-containing protein 81 (CCDC81)

<b>Product Code</b>	CSB-MP744397HU
<b>Abbreviation</b>	CCDC81
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ZN84
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLDTIARALQ DLGRQVLPTL PLSLQEEVSI IWGNVSEFVR RQLTLHKGVQ IPAFGTFTFI RQKLEVGNNK FILIQRPVFI MVEKLVQIHG LKQNKVYTPG EIIIVPLNFV MISLEGPVFN R DVVEGCVKET LLFLSRISM KQNVEFTFKG IGVLMIRDSK VKMRFYKDFL CTMDGSGALA KALANRPGTV DSVLSSREAL RKWPSSVLAF PRIELKEMEN KLPMETLVEE CGENRERKCK LKDQSDKEEG TRDISSPKRL RDRQALFPAK VTNVSLLEKF ERSESGGKIM TPESLSYPS LKHDSMMPKQ TSPACQDHNK AGQEMCYVCL QRAQRNSLLY YSEERRREIE DERLIQQYQM LKDQEALFRH QMKSLATREQ NQKNAAYNLG VAEAIRNHKN EKPEFYKSFL FDKRPLSPAL NALKQEEYSR SLLKQMDNRQ ENEIKQRQYR ELMDRLEQVQ LTEELAAQRA KFLKDKMEET QCYKRALDAQ IKNKPSRLPP FEPDSSEPIF GKNEGELMVE KQKREQNYMK HQLEAAANHK RKAILHQLVD QRRDLQMLQR TQREHLADRT AELERVNRVN QCLQEDWERS AAMKKQRDLE DKAFERASDK LFLLDQCEKY RRCKQCQRRT SNVGESNLWP LNKFLPGSRL LV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CCDC81
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 81
<b>Expression Region</b>	1-652
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.