



Recombinant Human Glutamate decarboxylase-like protein 1 (GADL1)

Product Code	CSB-YP744410HU
Abbreviation	GADL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZQY3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSDSDRQCP VDGIDIDQQEM IPSKKNVAVLV DGVVLNGPTT DAKAGEKFVE EACRLIMEEV VLKATDVNEK VCEWRPPEQL KQLLDLEM RD SGEPPHKLLE LCRDVIHYSV KTNHPRFFNQ LYAGLDYYSL VARFMTEALN PSVYTYEVSP VFLLVEEAVL KKMIEFIGWK EGDGIFNPGG SVSNMYAMNL ARYKYCPDIK EKGLSGSPRL ILFTSAECHY SMKKAASFLG IGTENVCFVE TDGRGKMIPE ELEKQVWQAR KEGAAPFLVC ATSGTTVLGA FDPDLDEIADI CERHSLWLHV DASWGG SLM SRKHRKLLHG IHRADSVAWN PHKMLMAGIQ CCALLVKDKS DLLKKCYSAK ASYLFQQDKF YDVSYDTGDK SIQCSRRPDA FKFWM TWKAL GTLGLEERVN RALALSRYLV DEIKKREGFK LLMEPEYANI CFWYIPPSLR EMEEGPEFWA KLNLVAPAIK ERM MKK GSLM LGYQPHRGKV NFFRQVVISP QVSREDMDFL LDEIDLLGKD M
Source	Yeast
Target Names	GADL1
Protein Names	Recommended name: Glutamate decarboxylase-like protein 1 EC= 4.1.1.-
Expression Region	1-521
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.