



Recombinant Human Centromere protein U (MLF1IP)

Product Code	CSB-EP744481HU-B
Abbreviation	MLF1IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q71F23
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPRGRRRPR PHRSEGARRS KNTLERTHSM KDKAGQKCKP IDVFDFPDNS DVSSIGRLGE NEKDEETYET FDPPLHSTAI YADEEEFSKH CGLSLSSTPP GKEAKRSSDT SGNEASEIES VKISAKKPGR KLRPISDDSE SIEESDTRRK VKSAEKISTQ RHEVIRTTAS SELSEKPAES VTSKKTGPLS AQPSVEKENL AIESQSKTQK KKGISHDKRK KSRSKAIGSD TSDIVHIWCP EGMKTSDIKE LNIVLPEFEK THLEHQQRIE SKVCKAAIAT FYVNVKEQFI KMLKESQMLT NLKRKNAMMI SDIEKKRQRM IEVQDELLRL EPQLKQLQTK YDELKERKSS LRNAAYFLSN LKQLYQDYSD VQAQEPNVKE TYDSSSLPAL LFKARTLLGA ESHLRNINHQ LEKLLDQG
Source	E.coli
Target Names	CENPU
Protein Names	Recommended name: Centromere protein U Short name= CENP-U Alternative name(s): Interphase centromere complex protein 24 KSHV latent nuclear antigen-interacting protein 1 MLF1-interacting protein Polo-box-interacting protein 1
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.