



Recombinant Bacillus cereus Uncharacterized RNA methyltransferase BCE_0363 (BCE_0363)

Product Code	CSB-BP744743BQN
Abbreviation	BCE_0363
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q73EJ5
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ATCC 10987 / NRS 248)
Purity	>85% (SDS-PAGE)
Sequence	MSTKMTPPVE KNEFIDVVFE DLTHDGAGVA KVKGYPIFVK NGLPGEEAQI KIIKVKKNFA FGRLMKLHTE SPYRKDAECP VYNQCGGCQL QHLT YEGQLQ AKEKQVRDVM QRIGGLSDVP VHPVLGMKNP WVYRNKAQVP IGEREGGLVA GFYRQGTHDI INMESCLIQA EENDTLIQEV KRICEKHGIS AYNEERNKGT LRHVMARYGQ VTGEIMLVFI TRTAELPNKK AIIEE IATKF PEVKSIVQNV NPKRTNVIFG DKTTVLYGSE YIYDFIGDIK FAISARSFYQ VNPEQTKVLY DKTLEYAKLD GNETVIDAYC GIGSISLFLA QKAKKVYGV EIVPEAIEDAK RNAALNNMTN AEFVGVEAEV VIPKWYKEGV IADTMVVDPP RKGCD EALLN TIIDMKPKRV VYVSCNPATL ARDLKVLEEG GYKTQEVQPV DMFPHTTHVE CVVLM SRVEK
Source	Baculovirus
Target Names	BCE_0363
Protein Names	Recommended name: Uncharacterized RNA methyltransferase BCE_0363 EC=2.1.1.-
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.