



# Recombinant Treponema denticola Glycerol kinase (glpK)

<b>Product Code</b>	CSB-EP744783TPP
<b>Abbreviation</b>	glpK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q73LE9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Treponema denticola (strain ATCC 35405 / CIP 103919 / DSM 14222)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKKYVLAFDQ GTTSCRILF DKDGKKIETA QQEFSQIFPK QGWVEHDAME IWGKQSGVAR EVLERSGVST QEIAAIGITN QRETTVVWNK NTGRPVCNAI VWQCRRTADI CDALKEKGLS DSIRKKTGLI IDAYFSGTKI KWILDNVPEA RTLAEKGELL FGNIDTWLIW NLTRGKVHVT DYTNASRLL FNIHTLQWDK ELLKAMDIPE SMLPEVKPSS YVYGYTEEHT FGGSKIPIAG AAGDQQAALF GQACFEEGSA KNTYGTGCFM LMNTGDKIIE SENGLTTIA FGIDNSVKYA LEGSSFIAGA AVQWLRDELK LIYTAHETAY YAGLVNDTNG VYFVPAFSGL GAPYWDMYAR GALLGLTRGA KREHIVRAVL EAIAYQTKDV LYAMEKDSKI NLKSLKVDGG ACANNFLMQF QSDILNVPVL RPYEKETTAL GAAYLAGLAV GFWKEQGEIK RIQDIEREF RPDMEEEKRKT LYAGWKKAVE RSMKWA
<b>Source</b>	E.coli
<b>Target Names</b>	glpK
<b>Protein Names</b>	Recommended name: Glycerol kinase EC= 2.7.1.30 Alternative name(s): ATP:glycerol 3-phosphotransferase Glycerokinase Short name= GK
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.