



Recombinant *Ashbya gossypii* Pre-mRNA-splicing factor SLT11 (SLT11)

Product Code	CSB-MP744951DOT
Abbreviation	SLT11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q750K9
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MDEETPSICE QCLTDADLRM TRAARGAECK ICTLPFTLYH FKPPGAPRVT KTLVCRRCAA QRNVCQCCML DLAWKLPVAL RDELVSLVQG SDERTPEASN EMVRRFLALR GGQLGGARLT ADSAPLRELM GKMRAVAAAA APAARAHGGR EDAGEVPETQ LPPGGALATP ASRSFFIYGV DPALPEWELV DAVSQLVRTP DWRDAGSVSV VVRQDARCAG IRFRREELAQ AFVARLETLP AAAGAPKGV L RIRHLRMHV V AWPEFHAAF GETKQQT EGI ARLATKTMRA DAEGASPARP AAPASAAGKV AKKAAAPRRQ GRRGPRPTLE L
Source	Mammalian cell
Target Names	SLT11
Protein Names	Recommended name: Pre-mRNA-splicing factor SLT11
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.