



Recombinant *Ashbya gossypii* pH-response transcription factor pacC/RIM101 (RIM101)

Product Code	CSB-EP744970DOT
Abbreviation	RIM101
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q753Y2
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSHINHLLNS DEAIPTALR SCAQRSPMSS DDELVGTAGY RENQPVPGST LPTSPSDVST DTENDEHQRL LCQWDACGSE FSQPELLYHH LCQDHVGRKS QRNLQLNCHW GSCTTKTVKR DHITSHLRVH VPLKPFSCST CSRKFKRPQD LKKHLKVHME DTMKERSRAA PGSRGVRKTG VNKGSALQEK ARTLPNLTVE SFVSEQEMQNY YPYYKSRQHL DETLSHILP PPAALGGTLA SEPPSYTRKA VSFFTTLSQD MSRRPLSLAP CNSPGPAGKM VMLPRPEQQY ARVPRYPAMP ELPPLVTSPG AESHALPRGH NFRPAPLIPG SMLPSLSRFS HDKSQFVARN SYSLNQKASL NEELCENELD IALENLALDD TDALRAELQT VNIKDYLTC LLEDVYTES FELPASVEKS HSTLRKYPQI KV
Source	E.coli
Target Names	RIM101
Protein Names	Recommended name: pH-response transcription factor pacC/RIM101
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.