



Recombinant *Ashbya gossypii* Ribosome biogenesis protein YTM1 (YTM1)

Product Code	CSB-EP744983DOT
Abbreviation	YTM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q756D0
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MVEDKSQVKV KFFTREQDES LHVLDAPLFA PVSLKRYGLS EIINHLLGLA EPVPFDLID GQLLRCSLQD YLTRKGLSSE AVLNVEYTRA VLPPSYLKSF SNEDWVSALD VGAERIVSGS YDGVVRTWNL SGKIEKQYSG HTGAVRAVKF ISSTRLVSGG NDRTRLRWKT KNDDVKHVDE LEGTEEAHTL AILEGHQAPV VSVDVQGDRI LSASYDNSIG FWSTNHKDMT AVDPMDSLGD KASSAAKKRR KLTMKDGSVR RRAPLSLLES HKAPVEQVIF ASNDSTVAYS VSQDHTIKTW DLVTSRCVDT KSTSYLLSM VELPKLRLLA CGSSARHITL HDPRADSSAK ITQQQLLGHK NFVVALDTCP ENEYMLCSAS HDGTVKVWDI RSSSSIYTIT REDQSVEKGI NDKVFAVKWA KGVGIISGGQ DKKIQFNKGD NIFKN
Source	E.coli
Target Names	YTM1
Protein Names	Recommended name: Ribosome biogenesis protein YTM1
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.