



Recombinant *Ashbya gossypii* rRNA 2'-O-methyltransferase fibrillar in (NOP1)

Product Code	CSB-EP744986DOT-B
Abbreviation	NOP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q756P0
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MAFQPGSRGG RGGARGGARG GARGGRGGFG GRGGSRGGRG GFDSRGGARG GFGGRGGSRG GPRGGPRGGA RGGRGGARGG AKGGAKVVIE PHKHAGVFIA RGKEDLLVTK NVAPGESVYG EKRI SVEEPA SEEGVPPTKV EYRVWNPFRS KLAAGIMGGL DELFIAPGKK VLYLGAASGT SVSHVADVVG PEGLVYAVEF SHRPGRELIS MAKKRPNVIP IEDARHPQK YRMLIGMVDA VFADVAQPDQ ARIIALNSHM FLKDQGGVVI SIKANCIDST VDAETVFARE VQKLREEKIK PLEQLTLEPY ERDHCIVIGR YMRSGL
Source	E.coli
Target Names	NOP1
Protein Names	Recommended name: rRNA 2'-O-methyltransferase fibrillar in EC= 2.1.1.-
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.