



Recombinant *Ashbya gossypii* Probable kynurenine formamidase (BNA7)

Product Code	CSB-BP745005DOT
Abbreviation	BNA7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q759T3
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MQRDYTATEV FQESGPDAIV YLHGGAWIDV RNSPRDFAEL AARVREDAPL ASQYAVTYRL SPEVQHPVHV ADCVEHISQL IVEKGIRKLH LLGHSVGATL CWQILTALPG DRRYPEGPDF AMKLALVRSS LAHVFLVDGI YSLRALLEEY PDYDYFVNKA FKSLKDFEDP SESAERISAG PILHVIHSYK DELLTLKQTQ YLTAVLTAHQ IPYKLYVDDL GLHEDVYRNE KVAAYIVAQL P
Source	Baculovirus
Target Names	BNA7
Protein Names	Recommended name: Probable kynurenine formamidase EC= 3.5.1.9 Alternative name(s): Biosynthesis of nicotinic acid protein 7 Probable N-formylkynurenine formamidase
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.