



# Recombinant *Ashbya gossypii* Glycylpeptide N-tetradecanoyltransferase (NMT1)

<b>Product Code</b>	CSB-YP745036DOT
<b>Abbreviation</b>	NMT1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q75EK2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGDSGDAQKM QRLLEMLALN SGDMSKLTQQ QRKAFEEYKF WKTQPVARFD EKVEEEGPIN PRRVEDVRD EPYPLLEEF WRTMDITTGQ DLEDVFNLLN ENYIEDKDST FRFNNTREFF NWALKPPGWR KEWHVGVVRV QSGRLVAFIS AVPTTLEVRG REMKSVEINF LCIHKKLSK RLAPILIKEI TRRVNKCDIW HALYSAGIVL PSPISTCRYT HRPLNWSKLF DVGFTALPAN ATKTQMLAKY TLPKKPLVEG LRPMTDADVD GAFDLFNRYQ KRFELIQTFD KSEFRHWFLG NEETPSVIYS YVQNSEGKI TDFVSFYSLP FTILKNPLHK ELGIGYLFYY ASDADFDYED RYDPTATELL RKRLTQLIND VCILARDLKM DVFNALTSQD NALFLEDLKF GPGDGFNLFY LFNYRANPIR GGLTEDKKFD AKNRSNQGCV ML
<b>Source</b>	Yeast
<b>Target Names</b>	NMT1
<b>Protein Names</b>	Recommended name: Glycylpeptide N-tetradecanoyltransferase EC= 2.3.1.97 Alternative name(s): Myristoyl-CoA:protein N-myristoyltransferase Short name=NMT Peptide N-myristoyltransferase
<b>Expression Region</b>	1-452
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.