



# Recombinant *Drosophila melanogaster* ATPase ASNA1 homolog (CG1598)

<b>Product Code</b>	CSB-EP745315DLU-B
<b>Abbreviation</b>	CG1598
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7JWD3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADNLEPLEP SLQNLVEQDS LKWIFVGGKG GVGKTTCCSSS LAVQLSKVRE SVLIISTDPA HNISDAFDQK FTKVPTKVNG FDNLFAMEID PNAGLNELPE EYFDGENEAL RVSKGVMQEM INALPGIDEA MSYAEVMKLV KGMNFSVVVF DTAPTGHTLR LIAFPQVVEK GLGKLLRLKM KVAPLLSQFV SMLGMADVNA DTLSQKLDDM LRVITQVNEQ FKNPDQTTFFV CVCIAEFFSL YETERLVQEL TKCGIDVHNI IVNQLLFLQN SHDSCSMCAS RFKIQEKYLD QIADLYEDFH VTKLPLLEKE VRGPESIRSF SENLMKPYNP KGEPKE
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	CG1598
<b>Protein Names</b>	Recommended name: ATPase ASNA1 homolog EC= 3.6.-.- Alternative name(s): Arsenical pump-driving ATPase homolog Arsenite-stimulated ATPase
<b>Expression Region</b>	1-336
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.