



Recombinant Drosophila melanogaster Proteasomal ubiquitin receptor ADRM1 homolog (CG13349)

Product Code	CSB-MP745318DLU
Abbreviation	CG13349
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7K2G1
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MFGRQSGGLGS SSNSSNLVEF RAGRMNMVGK MVHPDPRKGL VYMTQSDDGL MHFCWKDRTS GKVEDDLIVF PDDFEYKRVD QCKTGRVYVL KFKSSTRMF FWMQEPKTDK DDEQCRRINE LLNNPPSAHQ RGGGGSNDGD LQYMLNMSQ QQLMQLFGGV GQMGGSSLL GQMNSRTPSS RNTSSSGGGG ASALQTPENV SVPRTPSAPS KSGSSRSSN VNSQVGEAG SSVADAPGR SLNIDLSTAL PGADAINQII ADPEHVKTLL VHLPESEDVD DDRKQIQKDN ITSPQFQQAL AQFSSALQSA QLGPVIKQFE LSNEAVAAAF SGNLEDFVRA LEKSLPPGAT MGGKPSASEK KASDPETPTS VARDENTDPA TEKQEEKQK
Source	Mammalian cell
Target Names	Rpn13
Protein Names	Recommended name: Proteasomal ubiquitin receptor ADRM1 homolog Alternative name(s): p42E
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.