



Recombinant Flap endonuclease 1 (FEN1)

Product Code	CSB-BP745758PMT
Abbreviation	FEN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7RME3
Product Type	Recombinant Protein
Immunogen Species	Plasmodium yoelii yoelii
Purity	>85% (SDS-PAGE)
Sequence	<p> MGIKGLTKFI ADTAPNAIKE IKIENLMGRV VAIDASMSLY QFIIAIRDGD QYGNLMNESG ETTSHISGLM SRTIKLMENG LKPIYVFDGA PPELKGSELE KRGEKRQKAE ELLIKAKAEG NLEEIKKQSG RTVRVTKKQN EEAKLLTLM GIPVIESPCE AEAQCAFLTK YEMAHATATE DADALVFGTK ILIRNLNANA SSNKNKNKNS SKRGYILTEI NLEQVLKGLK LTMDEFIDFC ILCGCDYCDT IKGIGSKTAY NLIKEYNCIE NIIKNIDQNK YQVPANFKYV EARQSFINPK VLEKSEVKID WCEPKIEELK TFLIKEHNFN EVRVYNYITR LLKARKVTTQ RRLDTFFFTTC TTKSTKLIIE ESQKEVLKTK GKGKKREINN DSSTKLNAAK KKTNVKDEKK NNEKVDELKN KSDENLVKDE EDDQDDYDQN IFDEKTNSES GNIKNENVKE DVSSNDITMD IPKCTNDIVC </p>
Source	Baculovirus
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.