



Recombinant *Xenopus laevis* Inositol-tetrakisphosphate 1-kinase (itpk1)

Product Code	CSB-MP745821XBE
Abbreviation	itpk1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SY78
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MQTFLKGKRV GYWMSEKKIK KLNFAQAFADL CRKRGIEVVQ LNLAKPIEDQ GPLDVIIHKL TDVILEADQK DSESMQLVQR FQDYIEAHPE TIILDPLPAI RTLLDRSKSY ELIRRIETYM QDERICSPFF MELMAECDDED TLKILEKNGL AFPLVCKTRV AHGTNSHEMA IIFNPEGLWS IKPPCVIQSF ISHNAVLYKV FVVGESYTVV ERPSLKNFSL GASDRASIFF NSHNVSKPES SSVLTALEKV EGVFERPCDE VIRGISKALR QALGISLFGI DIIINNKTGQ HAVIDINAFP GYEGVPEFFT DLLNHITTIL QRPDQSANKS SVEQTGSSLG ERLCCTPPAR TDPWIVETDT SGSVKLQSQR LGCNSSVSPS FQQHCVASLA TKASSQ
Source	Mammalian cell
Target Names	itpk1
Protein Names	Recommended name: Inositol-tetrakisphosphate 1-kinase EC= 2.7.1.134 Alternative name(s): Inositol 1,3,4-trisphosphate 5/6-kinase Short name= Inositol-triphosphate 5/6-kinase Short name= Ins(1,3,4)P(3) 5/6-kinase EC= 2.7
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.