



Recombinant Acipenser baerii Keratin, type I cytoskeletal 18 (krt18)

Product Code	CSB-MP745822AAAJ
Abbreviation	krt18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SYF8
Product Type	Recombinant Protein
Immunogen Species	Acipenser baerii (Siberian sturgeon)
Purity	>85% (SDS-PAGE)
Sequence	SYRPGSYSV SSMRPVGSVR SSQVMTVQRS MPLASAASVY GGAGGRGSRI SVGGSSSGFG SGLGSGAGGS YSSMSVSGSG LVGNEKETMI GLNDRLAAYL ETVRNLEQAN SKLEFQIREA LEKKGPTTRD LSPFEKTLED LRKKVYDMTM DNSRLVLQID NSRLAADDFR VKFESEYSIR QSVESDIIGL RKVIDDTNMG RMNLESEIES LKEELIFIKQ NHNQEVNDLR NQIAQSGVQV DVDAPKGQDL AQVLAEVRAQ YESMAQKNRD ELKAWHENKL TEVEVEVIQN TEALQGARTE VTELRRQM QS LEIELESQRS MKASLEDSL R DTEMRNNMEM ERYNNMILQL EAELGQLRGN IQMQASEYEA LLNIKMKLEA EIATYRLLD GEDFRLQDAL VDQSSTKSIK KVTVTQTLVD GKVVSESTNT KEIGK
Source	Mammalian cell
Target Names	krt18
Protein Names	Recommended name: Keratin, type I cytoskeletal 18 Alternative name(s): Cytokeratin-18 Short name= CK-18 Keratin type Is Keratin-18 Short name= K18
Expression Region	2-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.