



Recombinant Danio rerio Cytosolic 5'-nucleotidase 3 (nt5c3)

Product Code	CSB-YP745823DIL
Abbreviation	nt5c3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SYN4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MMPEFEKNTV HIRDPERVEQ IICGLIKGGA SKLQIITDFD MTLSRFAVNG KRCPSCHNII DNSKLVTDCC RKKLVHLKET YYPIDPHL TMEEKYPFMV EWYFKSHTLL VEQRLEKDKL PEAVRESVVS LKEGYEQFFD RLHQHSVPVF IFSAGLGDVL EEIRQAGVY HPNVKVVS NF MDFDDNGVLK GFKGELIHVY NKHDGALRNT EYFKQLKDNG NIILLGDSL G DLT MADGVPN VENILKIGYL NDKVEELLEK YMDSYNIVLA RDETLEVPNS ILQKIL
Source	Yeast
Target Names	nt5c3
Protein Names	Recommended name: Cytosolic 5'-nucleotidase 3 EC= 3.1.3.5 Alternative name(s): Cytosolic 5'-nucleotidase III Short name= cN-III
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.