



Recombinant *Xenopus tropicalis* Eukaryotic translation initiation factor 2 subunit 1 (eif2s1)

Product Code	CSB-EP746783XBF-B
Abbreviation	eif2s1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GL89
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	PALSCRFYQ HKFPEVDDVV MVNVRSAEM GAYVSLLEYN NIEGMILLSE LSRRRIRSIN KLIRIGRNEC VVVIRVDKDK GYIDLSKRRV SPEALKCED KFTKSKTVYS ILRHVAEVLE YTKDEQLESL YQRTAWVFDE KYKKPGYGAY DAFKHAVSDP DVLDGLDLNE DERRVLIDNI NRRLTPQAVK IRADIEVACY GYEGIDAVKD ALRAGLSCST ENMPIKINLI APPRYVMTTTLERTEGLSV LNQAMSVIKE RIEEKRGVFN VQMEPKVVTD TDETELARQL ERLEKENAEV DGDDDAEEME AKTED
Source	E.coli
Target Names	eif2s1
Protein Names	Recommended name: Eukaryotic translation initiation factor 2 subunit 1 Alternative name(s): Eukaryotic translation initiation factor 2 subunit alpha Short name= eIF-2-alpha Short name= eIF-2A Short name= eIF-2alpha
Expression Region	2-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.