



Recombinant Human THO complex subunit 7 homolog (THOC7)

Product Code	CSB-MP746911HU
Abbreviation	THOC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6I9Y2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GAVTDDEVI RKRLIDGDG AGDDRRINLL VKSFIKWCNS GSQEEGY SQY QRMLSTLSQC EFSMGKTLV YDMNLREMEN YEKIYKEIEC SIAGAHEKIA ECKKQILQAK RIRKNRQEYD ALAKVIQHHP DRHETLKELE ALGKELEHLS HIKESVEDKL ELRRKQFHVL LSTIHELQQT LENDEKLSEV EEAQEASMET DPKP
Source	Mammalian cell
Target Names	THOC7
Protein Names	Recommended name: THO complex subunit 7 homolog Alternative name(s): Functional spliceosome-associated protein 24 Short name= fSAP24 Ngg1-interacting factor 3-like protein 1-binding protein 1 Short name= NIF3L1-binding protein 1
Expression Region	2-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.