



Recombinant Danio rerio SUMO-activating enzyme subunit 1 (sae1)

Product Code	CSB-BP746933DIL
Abbreviation	sae1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6IQS6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MIDTIEKEDT IISEEEAAQY DRQIRLWGLD AQKRLRGSRV LLVGLRGLGA EVAKNLILAG VKGLTLLDHE QVTEESRRAQ FLIPVDADGQ NHAQASLERA QFLNPMVEVK ADTEPVESKP DDFFFQFQDAV CLTRCSRDLM VRVDQLCASR NIKVFCGDVY GYNGYMFSDL GQEYHYVEEK PKVVKGSNEA NDGPEAKKPK IDPNETTMVK KTISFCSLKE ALEVDWTEK AKSSLKRIPA DYFLLQVLLK FRTDKGRDPQ PDSFAEDSQL LLQIRDDVLE TMGLSSDLLP NTFVSYCFSE MSPVCAVVG VLGQEIVKAL SQRDAPHRNF FFFDGLKGS VVDYFSSK
Source	Baculovirus
Target Names	sae1
Protein Names	Recommended name: SUMO-activating enzyme subunit 1 Alternative name(s): Ubiquitin-like 1-activating enzyme E1A
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.