



Recombinant Danio rerio Homeobox protein Hox-C13b (hoxc13b)

Product Code	CSB-EP746954DIL-B
Abbreviation	hoxc13b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6JIY4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEGLSGNCPA SHCRDFISHP ALGRHSGSLA SHQGTVYPDI TTQDAGRQFP APQASSGTSL GYGYAFGSPY YGCRLSYSHN VNLQQKSYHP AEKYMETS GA LPAEELSSRS KEFAIYPSFA SSYQTVPGYL DVPVVPGIS A HPESRHEALF PMDSYQHWAL SNGWDEQLYC SKEQTHFNHL WKSQFSDVVP HQAEMNGYRR GRKKRVPYTK IQLKELEKEY AASKFITKDR RRRISATTSL SERQVTIWFQ NRRVKEKKFV CKSSKSSTNM HSV
Source	E.coli
Target Names	hoxc13b
Protein Names	Recommended name: Homeobox protein Hox-C13b
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.