



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-F2-2 (CYCF2-2)

Product Code	CSB-BP746963OFG
Abbreviation	CYCF2-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6K8S0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDSIMEPYVA DLLADDITAS MVELLPGDGG AAQMDVGVLD AYLRAIGALP AHPAAPGADL AAAAEVESMA SNDDTNGVLY DWDTKVDVKV PCALLPPPPG FPPLPVPGLA DEPVYAAPAR HLPPPPGFPP LPVPGLADEP VYAAPARRLP PPPGFPPLPV PAKAEPVYAA PVDEGDaira FMQQLEWSEQ YNGDNDAPAP DNSTASRPQL CAPYDDDIDA NLRDMEKDAA QRPSPDYLDT VQGGQISAAA RASLVAWMGR LTHRYELAAG TLHRAVSyFD RFLSVRALPS YTAHQLSLVA ATAVYTAAY EDQGTvFKLD AREIASYGEF ASAQEVLAME REMMAALGYR LGGPNAETfV EHFTRYSKGK EELRVQRLAC HVADRSLESY GCLGYLPSMV AAA AISIARW TLNPPGALPW SSELQELTGY SSQDISSCIL TVLNTQ
Source	Baculovirus
Target Names	CYCF2-2
Protein Names	Recommended name: Cyclin-F2-2 Short name= CycF2;2
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.