



Recombinant Arabidopsis thaliana AMSH-like ubiquitin thioesterase 2 (AMSH2)

Product Code	CSB-EP747311DOA
Abbreviation	AMSH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NKP9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVTLLSSPSPS LSCVENVTCK SSHVSRVLIS GTDNINHGES SEAKILRDVH ISERLLEDFT ELARENTEKD LETCGTAAAF LERGIFYVTT LIIPKQUESTS NSCQAMNEVE VFSIQNEREL YPVGWIHHP SQGCFMSSVD LTHYSYQVM VPEAFAIVVA PTDSSKSYGI FKLTDPGGME VLRGCSETGF HPHKEPEDGN PVYEHCSNVY KNSNLRFEIF DLR
Source	E.coli
Target Names	AMSH2
Protein Names	Recommended name: AMSH-like ubiquitin thioesterase 2 EC= 3.4.19.- Alternative name(s): Deubiquitinating enzyme AMSH2
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.