



Recombinant Human 3-hydroxyisobutyryl-CoA hydrolase, mitochondrial (HIBCH)

| | |
|--------------------------|--|
| Product Code | CSB-YP747360HU |
| Abbreviation | HIBCH |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6NVY1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | DAAEEVLL EKKGCTGVIT LNRPKFLNAL TLNMIRQIYP QLKKWEQDPE TFLIIKGAG GKAFCAAGDI RWISEAEKAK QKIAPVFFRE EYMLNNAVGS CQKPYVALIH GITMGGGVGL SVHGQFRVAT EKCLFAMPET AIGLFPDVGG GYFLPRLQ GK LGYFLALTGF RLKGRDVYRA GIATHFVDSE KLAMLEEDLL ALKSPSKENI ASVLENYHTE SKIDRDKSFI LEEHMDKINS CFSANTVEEI IENLQQDGSS FALEQLKVIN KMSPTSLKIT LRQLMEGSSK TLQEVLTMEY RLSQACMRGH DFHEGVRAVL IDKDQSPKWK PADLKEVTEE DLNNHFKSLG SSDLKF |
| Source | Yeast |
| Target Names | HIBCH |
| Protein Names | Recommended name: 3-hydroxyisobutyryl-CoA hydrolase, mitochondrial EC=3.1.2.4 Alternative name(s): 3-hydroxyisobutyryl-coenzyme A hydrolase Short name= HIB-CoA hydrolase Short name= HIBYL-CoA-H |
| Expression Region | 33-386 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.