



Recombinant Human Armadillo repeat-containing protein 6 (ARMC6)

Product Code	CSB-MP747367HU
Abbreviation	ARMC6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NXE6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSERCCSRYS SGASIGCTPT STQAKMVSKR IAQETFDAAV RENIEEFAMG PEEAVKEAVE QFESQGV DLS NIVKTAPKVS ADGSQEPTH D ILQMLSDLQE SVASSRPQEV SAYLTRFCDQ CKQDKACRFL AAQKGAYPII FTAWKLATAG DQGLLLQSLN ALSVLTGQ P DLLDAQGLQL LVATLTQNAD EADLTCSGIR CVRHACLKHE QNRQDLVKAG VLPLLTGAIT HHGHHTDVVR EACWALRVMT FDDDIRVPFG HAHNHAKMIV QENKGLK VLI EATKAFLDNP GILSELGTL SRLAIRNEFC QEVVDLGGLS ILVSLLADCN DHQMRDQSGV QELVKQVLST LRAIAGNDDV KDAIVRAGGT ESIVAAMTQH LTSPQVCEQS CAALCFLALR KPDNSRIIVE GGGAVAALQA MKAHPQKAGV QKQACMLIRN LVAHGQAFSK PILDLGAEAL IMQARSAHRD CEDVAKAALR DLGCHVELRE LWTGQQRGNLA P
Source	Mammalian cell
Target Names	ARMC6
Protein Names	Recommended name: Armadillo repeat-containing protein 6
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The function of this gene s protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.